


Red Light Therapy Wrap

Instructional Manual

Red-Light-Therapy for Human & Animals !

Thank you for purchasing our red light , we hope you can reap benefits of red light therapy at home .



Say hello to your flexible red light wrap !

Before starting your first session , please make sure to read this manual as it may explain important operating procedures .

Notes : Each led diode inside with three chips at 1 pc 660nm and 2 pcs 850nm .

This package includes

- 1* Red light therapy pad
- 1* Power cord
- 1* Adaptor
- 1* Elastic Strap
- 1* User Manual

Specifications

| | |
|--------------|-----------------------------------|
| Rating Power | 175W |
| Actual Power | 45W |
| Product Size | 890*280*6.3MM |
| LED QTY | 350*0.5W |
| Net Weight | 0.44 KG |
| Wavelength | 1*660nm ; 2*850nm |
| EMF | 0 UT |
| Irradiance | 0 Inch: 124.5 mW/c m ² |
| Timer | 20 Minutes |
| Material | Fabric |
| Lifespan | Over 50,000 Hours |

The red light device uses combo red 660nm and Infrared 850nm leds which is chosen to give the ideal wavelengths found in the scientific research of the last three decades .

660nm of red light

This light can be seen by the human eyes and targets skin cells

850nm of near infrared light

This light cannot be seen by the human eyes and targets deep tissue

User Instruction

- 1) Plug into the AC socket and the button illuminates green than press the power button to turn the device on.
- 2) Distance: 0 inches for best results, place anywhere on the body you seek relief.
- 3) Select options:

Press the power button once quickly so the power button illuminates green, then:

- a. Press the button once, and all LEDs turn on.
- b. Press again, and only the red LEDs will be on, and near infrared off.
- c. Press again, only near infrared LED on, and red LEDs off. **Note that infrared light is invisible to the human eye so it may appear that the LEDs are off when in fact the infrared is on.**
- d. Press again, and all LED diodes will turn off.



If you press the power button for more than 3 seconds, the power button illuminates blue, and:

- a. Press the button once, and all LEDs come to pulse mode.
- b. Press again, and only the red LEDs come to pulse mode.
- c. Press again, and only the near infrared LEDs come to pulse mode.
- d. Press again, and all LED diodes will turn off.



4) Usage Frequency: Use device 1-2 times per day for at least 4 weeks, and for a duration of 20 minutes on one treatment area before moving the device to the next treatment area.

5) Unplug device and store in cool dry area.

6) Maintenance: For cleaning we recommend using a damp towel or toothbrush to **slightly** clean your wrap and elastic strap. Make sure the device is **unplugged** from any power, cords disconnected and turned off. Only clean the fabric portion of the wrap and not the LED diodes directly. Do not submerge the wrap in water. Allow the wrap to air dry and make sure to keep it out of direct sunlight.

Benefits backed by science :

Research shows , that the wavelengths in the spectrum of 660-900nm show optimal biotical response .

- **Muscle and Joints pain**
- **Speed wound healing**
- **Hair Re-growth**
- **Enhanced blood circulation**
- **Aging spots and wrinkle removal**
- **Increase collagen**
- **Weight loss**
- **Anti-inflammatory**
- **Decreases pain**

Even though this technology don't have side effects , we suggest taking the following precautions :

- ▷ **This device is water-resistant but not waterproof ! Do not submerge the wrap or its components in water.**
- ▷ **Not intended for the use by children and pregnant woman.**
- ▷ **Do not exceed 30 Min. on one same treatment area.**
- ▷ **Don't stare directly at the light .**
- ▷ **For cleaning instruction please wipe the wrap slightly with wet tissue or dry towel .**

Application

Here are some examples of how you can use the red light device . You and your pets can enjoy the red light at home .



