


WIRELESS Rectal Repro B/W Ultrasound X Probes



X Probe WIRELESS Rectal Repro B/W Ultrasound Scanner



X Probe-L – Linear Probe
with a radius of 64 mm with 7.5 or 10.0 MHz
and a Scan Depth of 100 mm
Suited for small placenta in cows



X Probe-C – Convex Probe
with a radius of 60 mm with 3.5 or 5.0 MHz
and a Scan Depth of 200 mm
Suited for larger placentas in cows



X Probe-N – Linear Probe
with a Probe length of 20 mm with 7.5 or 10.0 MHz
and a Scan Depth of 100 mm
Suited for ewes and nanny's

Specifications:

- WIRELESS connectivity from the wand to your iOS (Apple) device
- Working System: iOS (Apple)
- Battery Working Time: 6 to 8 hours
- Image save forma: jpg, avi, and DICOM
- Image Adjust: Gain, DYN, Focus, Depth, Harmonic, & Denoise
- Image Mode: B, B/M
- Measure: Length, Area, Angle, & Obstetrics
- Cineplay: Auto or Manual, Frames can be set as 100/200/500/1000
- Image Frame Rate: 18 frames per second
- Average Weight – 1 ½ lbs.
- Charge: by USB

