

NAVIGATOR AI

Stationary X-Ray System

Compatible with any flat-panel detector & acquisition software, it's a perfect blend of advanced technology, user-friendly design, and compatibility, making it an ideal solution for efficient and accurate veterinary diagnostics.

**5 YEAR
WARRANTY**



Streamlined Workflow



Stand-Alone Operation



User-Friendly Interface



Universal Compatibility



Dynamic Exposure Adjustment



Patented Automatic Thickness Measurement



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983


www.vetxray.com

Contact us:

605.341.2433
scott@vetxray.com
2325 E. Saint Charles St.
Rapid City, SD 57703



NAVIGATOR AI



Key Features:



Distinctive Panda Design



Universal Compatibility



Four-Way Floating Tabletop



Restraining Straps



Urine Track



Auto-Tube Calibration



Canon E7239X Tube



Max Amperage input: 16 amps



32kw, single phase high frequency generator



Patented Automatic Thickness Measurement



Collimator Field Illumination During Tabletop Movement



Easy Transportation, Mobility, and Installation



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com

Contact us:

605.341.2433

scott@vetxray.com

2325 E. Saint Charles St.
Rapid City, SD 57703

