

NAVIGATORX

Stationary X-Ray System

Analog X-ray systems provide a cost-effective solution for pet clinics, emphasizing functionality without the digital price.

**5 YEAR
WARRANTY**

- 32kW/400mA Compact HF Generator
- Hands Free Collimator light control
- Variable SID for table top and cassette tray imaging
- Rotating Tube Head and collimator
- X-ray tube Rotor Brake
- 4 Way Float Top Table – 5ft Standard
- 208/230/240 VAC Single Phase



ABOUT Diagnostic Imaging Systems

Diagnostic Imaging Systems (DIS) has been providing quality imaging products at the best possible possible prices since 1983. DIS products are specifically tailored to match all your practice needs.



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983



www.vetxray.com

Contact us:

605.341.2433
scott@vetxray.com
2325 E. Saint Charles St.
Rapid City, SD 57703



NAVIGATORX



Standard X-Ray Tube:

- 1.0/2.0 Focal Spots
- 16° Target angle
- 140,000 HU

Grid Cabinet:

- 8:1 103 line pair grid
- Stainless steel cassette tray

Options:

- Custom (38" to 78") length top
- Fractional Focal Spot X-ray tubes
- 3 phase and battery powered generators



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com

Contact us:

605.341.2433
scott@vetxray.com
2325 E. Saint Charles St.
Rapid City, SD 57703

