


Companion Animal X-Ray Package

Complete X-Ray System


Cost-effective solution for companion animal clinics, complete with stationary x-ray table and 17x17 cassette-size flat panel detector.



- 17x17 Cassette-size Flat Panel Detector
- Advanced Image Processing with Voyance Software
- 32kW/400mA Compact HF Generator
- Hands Free Collimator light control
- Variable SID for table top and cassette tray imaging
- Rotating Tube Head and collimator
- X-ray tube Rotor Brake
- 4 Way Float Top Table – 5ft Standard
- 208/230/240 VAC Single Phase

ABOUT Diagnostic Imaging Systems

Diagnostic Imaging Systems (DIS) has been providing quality imaging products at the best possible possible prices since 1983. DIS products are specifically tailored to match all your practice needs.



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com


Contact us:

605.341.2433
scott@vetxray.com
2325 E. Saint Charles St.
Rapid City, SD 57703



DR ^{17x17} wizard X

17"x17" Tethered Cassette-Size DR Flat Panel System



Advanced Voyance Software

- Superior Image Quality
- Exam-specific algorithms image processing
- Automatic image optimization
- Edit existing patient information/search patient by patient name, patient ID, Accession number, bar code, etc.
- True size for printing
- Resend / Reprint image
- Send/Print queue management
- User-friendly touch GUI
- Useful functions: Auto Cropping, Auto Grid Suppression, Auto Labeling, Auto Rotation, Free Rotation, Multi Layout, etc.
- Fully DICOM 3.0 Compatible: MWL, Send, Print, MPPS, etc.
- Customized functions and various themes

Specifications

Conversion Screen (Scintillator)	CsI / Gadox
Detector Type	Amorphous Silicon
Pixel Pitch	140 µm
Image Data (AD conversion)	16 bit
Image Resolution	3,072 x 3,072
Image Size	17" x 17"
Data Interface	Wired
X-Ray Sync	Lossless AED/AWC, Sync Trigger

Standard Vet Analog X-Ray System

Standard X-Ray Tube:

- 1.0/2.0 Focal Spots
- 16° Target angle
- 140,000 HU

Grid Cabinet:

- 8:1 103 line pair grid
- Stainless steel cassette tray

Options:

- Custom (38" to 78") length top
- Fractional Focal Spot X-ray tubes
- 3 phase and battery powered generators



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com

Contact us:

605.341.2433
scott@vetxray.com
2325 E. Saint Charles St.
Rapid City, SD 57703

