

# ELITE 125100

## Mixed Animal Portable X-Ray Unit

The Elite 125100 is an all-in-one portable X-ray machine that includes a 5 kW X-ray power unit. With a compact size, light weight, and high power density, it is ideally suited for a variety of portable X-ray applications in the clinic or in the field.



# 5kW




### ELITE Design

Compact and lightweight design for portable use.



### ELITE Connectivity

Easily connect to tablet through Bluetooth and adjust x-ray-parameter in real time.



### ELITE Accuracy

6 sets of APR exposure parameter shortcut keys to reduce errors and improve efficiency.

### ABOUT

### Diagnostic Imaging Systems

Diagnostic Imaging Systems (DIS) has been providing quality imaging products at the best possible prices since 1983. DIS products are specifically tailored to match all your practice needs.



# DIAGNOSTIC IMAGING SYSTEMS

QUALITY | INTEGRITY | VALUE...privately owned since 1983

[www.vetxray.com](http://www.vetxray.com)

Contact us:


605.341.2433

[scott@vetxray.com](mailto:scott@vetxray.com)

2325 E. Saint Charles St.  
Rapid City, SD 57703



# ELITE 125100



## Specifications

<b>Max Output Power</b>		<b>5 kW</b>
<b>Power Supply</b>		<b>AC 220V / 110V ± 10% 50 / 60 Hz</b>
<b>kV Range</b>		<b>40-125</b>
<b>MA's Range</b>		<b>1-100mAs</b>
<b>mA</b>		<b>10 - 100</b>
<b>Communication Mode</b>		<b>RS232/Bluetooth</b>
<b>SID Measurement</b>		<b>Ultrasonic</b>
<b>Light Source On Time</b>		<b>30 seconds</b>
<b>X ray Tube</b>	<b>Focal Spot</b>	<b>.6/1.8mm</b>
	<b>Target Angle</b>	<b>15 degree</b>
	<b>Small Focal Spot Power (kW)</b>	<b>1.27</b>
	<b>Large Focal Spot Power (kW)</b>	<b>5.3</b>
	<b>Inherent Filtration (mm Al @ 75 kV)</b>	<b>Min.0.65</b>
	<b>Anode Heat Storage (kJ)</b>	<b>30</b>
<b>Weight</b>		<b>35 lbs</b>



**DIAGNOSTIC  
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

[www.vetxray.com](http://www.vetxray.com)

Contact us:

605.341.2433  
scott@vetxray.com  
2325 E. Saint Charles St.  
Rapid City, SD 57703

