

ELITE Dental

Handheld Intra-oral X-ray System

The ELITE Dental has 7-inch LCD touch screen and Android Operating system. Images are transferred directly via USB - Realize, Shoot and See.



- Shoot and See: image enhancement, move, reverse color, contrast etc.
- Android Operating system: Easy to understand
- 7-inch LCD touch screen
- Independent patient management
- Automatic Pitch angle display
- Replaceable battery
- Light weight
- Support multi-language



**Images appear directly
on generator screen**



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983


www.vetxray.com

Contact us:

605.341.2433
scott@vetxray.com
2325 E. Saint Charles St.
Rapid City, SD 57703



ELITE Dental



Specifications

Operating System	Android
Screen	7-inch Touch Screen
kV	70
mA	2
Exposure Time (s)	0.04 - 2
Tube Model	0.4
Anode Angle	12°
Focal Spot	0.4mm
Data Transmission	USB
Total Filter	≥1.5mm AL
Weight	4.4 lbs
Internal Radiation Shielding	0.1Gy/h@ 1m
Operating Temperature	10 ~ 40
Storage & Transport Temperature with Package	-10 ~ 45
Operating Humidity (% RH) (Non-Condensing)	30 ~ 75
Storage & Transport Humidity with Package (% RH) (Non-Condensing)	30 ~ 75



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com

Contact us:

605.341.2433
 scott@vetxray.com
 2325 E. Saint Charles St.
 Rapid City, SD 57703

