

ELITE Size 4

Dental & Small Parts Sensor

Meet the ALL NEW Elite Size 4 Dental and Small Parts Sensor – perfect for clear, detailed imaging with advanced APS CMOS tech. Slim, CE-certified, and designed for comfort.



- Size 4 Sensor
- USB 2.0 Connection
- Enhanced Versatility
- Fast, High Quality Images
- More comprehensive area of imaging
- Efficient diagnoses of a wider range of patients



DIAGNOSTIC
IMAGING SYSTEMS

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com

Contact us:


605.341.2433

scott@vetxray.com

2325 E. Saint Charles St.
Rapid City, SD 57703



ELITE Size 4 Dental & Small Parts Sensor



Specifications

Quality Certification	CE
Instrument classification	Class II
Product Specification	Size 4
Sensor type	APS CMOS
External dimensions (mm)	77.1 x 53.8
Effective area (mm)	67.3 x 46.7
Sensor thickness	6.5mm
Gray scale range	0~16383
Power Supply	5.0 V (≥4.25 V)
Interface	USB2.0



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com

Contact us:

605.341.2433
 scott@vetxray.com
 2325 E. Saint Charles St.
 Rapid City, SD 57703

