

DR **Equine Wizard** PRO

10"x13" Wireless Equine Csl DR Flat Panel System

The DR Wizard Equine Pro is a versatile portable DR system that is easy to transport in the self powered, portable case without worrying about power cords to operate.



Lightweight High-Tech Composite

Extremely durable, lightweight portable veterinary DR case made with the latest high-tech composite material with shock absorbing edge reinforcement. The case is also splash proof (IPX4 rated).



Built-In X-Ray PC Workstation

Built-in high performance PC workstation with Windows operating system and AccuVue Vet software. Also comes with a built-in keyboard and trackpad mouse.



Padded DR Compartment

Easy removal and storage of the x-ray flat panel detector. Fits either panel size and protection case.

Huge 22" HD Touchscreen Monitor.

- Features a large, anti-glare tilting 22" touchscreen display when taking images on the go as well as a built in trackpad mouse for precise post processing tools.
- The integrated monitor is brought into the working position with the flick of your wrist
- Utilizes Power Saving Mode - Fold down the monitor to activate power saving mode and continue where you left off within a set timeframe. System auto shuts down if not reopened.
- Intuitive AccuVue Vet x-ray imaging software



www.vetxray.com

Contact us:

605.341.2433
scott@vetxray.com
2325 E. Saint Charles St.
Rapid City, SD 57703



DIAGNOSTIC IMAGING SYSTEMS

QUALITY | INTEGRITY | VALUE...privately owned since 1983



DR *Equine* Wizard PRO

Protection Cases

Available with both gridded and non-gridded protection cases to provide ultimate panel protection. Cases are equipped with an ergonomic handle and locking latch for secure handling.



Advanced Accu-Vet Software

Powerful engine offers a convenient user-friendly interface to display various anatomical tuning features from a single exposure to offer the most out of a diagnostic imaging solution.



Padded DR Compartment

Easy removal and storage of the x-ray detector. Fits panel and protection case.

Power Saving Mode

Fold down the monitor to activate power saving mode and continue where you left off with a set timeframe.



Intelligent Battery Design

unique battery docking system that allows you to swap in power tool rechargeable batteries on the go, to keep the case ready for clinical use for hours without recharging.



Specifications

Extremely durable case with edge shock absorption, splash-proof (IPX4), and is lockable	Comes with AccuVue Vet Software
Anti-glare 22" full HD touchscreen monitor, keyboard, and trackpad mouse	Includes DR protection case
Carrying strap for comfortable transport	(2) 20V/5 Ah Dewalt rechargeable batteries
Padded DR compartment for 10x13 detector	Use with your smartphone for wireless control
Connections: 1 x USB, 1 x LAN, 1 x WLAN, 1 x Bluetooth, 1 x DC socket	High performance PC with Windows OS
Total weight: 21 lbs (includes complete x-ray case, monitor, and detector)	Dimension: 21" W x 24" H x 20" D



Intuitive Operation with AccuVue Vet Software

- Modern graphical user interface (GUI) available in almost any language, with touchscreen operation to ensure a quick, efficient, and smooth workflow
- Capture of patient data via DICOM Worklist, BDT GDT HL / 7 or other protocols - data may also be captured manually
- Use of DICOM Procedure Codes for the transfer of all relevant examination data directly from the connected patient management system
- Freely configurable body parts with more than 200 projections and numerous possible adjustments already included
- Safe and fast registration of emergency patients allowing user to switch between projections to avoid having to re-position the patient frequently
- Allows the user to subsequently add images to an examination, even after that examination has been completed



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com

Contact us:

605.341.2433
scott@vetxray.com
2325 E. Saint Charles St.
Rapid City, SD 57703

