

ELITE 2kW-BT

Equine Battery-Powered Portable X-Ray Unit

The Elite 2kW-BT is built specifically for equine veterinarians, offering high-quality imaging in a lightweight, compact design. Its long-lasting battery and splash-proof construction ensure reliable performance in the field.



ELITE Portability

Lightweight, compact, and splash-proof design for easy portability in the field.



ELITE Battery Life

Long-lasting battery allows over 200 exposures at full power—ideal for all-day equine imaging.



ELITE Connectivity

App-based exposure control via Wi-Fi for convenient, remote operation.

ABOUT Diagnostic Imaging Systems

Diagnostic Imaging Systems (DIS) has been providing quality imaging products at the best possible prices since 1983. DIS products are specifically tailored to match all your practice needs.



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com

Contact us:

605.341.2433

scott@vetxray.com

2325 E. Saint Charles St.
Rapid City, SD 57703



ELITE 2kW-BT

Eco-friendly and User-oriented

- Multi-level leak-proof design ensures high image quality with low power and low radiation intensity.
- Compact and splash proof design allows operation in different environment
- APP exposure control through WIFI



Battery Power Supply

- A full charge enables over 200 times exposure at full power and over 400 times at average power.
- Normal exposure is available even when the power is lower than 20%.




Detachable Battery

Convenient Maintenance & Updates

- Exclusive USB update function supports instant system update
- Different prompt codes indicate system status



5 modes of exposure



App




Hand Switch*



Remote Control



Foot Switch



Delay Key

*Hand Switch is Standard, Others are optional accessories.

Specifications

Max Output Power	2 kW
kV Range	40-100kV
MA's Range	0.4-50mAs
mA	20/25mA



DIAGNOSTIC
IMAGING SYSTEMS

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com

Contact us:

605.341.2433
scott@vetxray.com
2325 E. Saint Charles St.
Rapid City, SD 57703

