

ELITE 3.2kW-BT


Mixed Animal Battery-Powered Portable X-Ray Unit

The Elite 3.2kW-BT is ideal for mixed animal practices, combining portability with enhanced power for versatile imaging across species. Its user-friendly design and app-based exposure control streamline diagnostics in any setting.




ELITE Balance

Balanced power and portability for versatile imaging across small and large animals.



ELITE Accuracy

High-quality images with low radiation exposure.



ELITE Maintenance

USB system updates and diagnostic codes simplify maintenance and system management.

ABOUT

Diagnostic Imaging Systems

Diagnostic Imaging Systems (DIS) has been providing quality imaging products at the best possible prices since 1983. DIS products are specifically tailored to match all your practice needs.



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com

Contact us:

605.341.2433

scott@vetxray.com

2325 E. Saint Charles St.
Rapid City, SD 57703



ELITE 3.2kW-BT

Eco-friendly and User-oriented

- Multi-level leak-proof design ensures high image quality with low power and low radiation intensity.
- Compact and splash proof design allows operation in different environment
- APP exposure control through WIFI



Battery Power Supply

- A full charge enables over 200 times exposure at full power and over 400 times at average power.
- Normal exposure is available even when the power is lower than 20%.




Detachable Battery

Convenient Maintenance & Updates

- Exclusive USB update function supports instant system update
- Different prompt codes indicate system status



5 modes of exposure



App



Hand Switch*



Remote Control



Foot Switch



Delay Key

*Hand Switch is Standard, Others are optional accessories.

Specifications

Max Output Power	3.2 kW
kV Range	40-125kV
MAs Range	0.4-200mAs
mA	10-80mA

