

ELITE 5kW-BT

Mixed Animal Battery-Powered Portable X-Ray Unit

The Elite 5kW-BT delivers maximum power and performance for high-volume, mixed animal practices. Engineered for durability and efficiency, it provides superior image quality with extended battery life for demanding clinical workflows.



ELITE Power

Maximum power output for rapid, high-resolution imaging in busy veterinary settings.



ELITE Durability

Durable, splash-proof build designed for use in diverse and demanding environments.



ELITE Battery Life

Extended battery life delivers consistent performance, even below 20% charge.

ABOUT

Diagnostic Imaging Systems

Diagnostic Imaging Systems (DIS) has been providing quality imaging products at the best possible possible prices since 1983. DIS products are specifically tailored to match all your practice needs.



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com

Contact us:

605.341.2433

scott@vetxray.com

2325 E. Saint Charles St.
Rapid City, SD 57703



ELITE 5kW-BT

Eco-friendly and User-oriented

- Multi-level leak-proof design ensures high image quality with low power and low radiation intensity.
- Compact and splash proof design allows operation in different environment
- APP exposure control through WIFI



Battery Power Supply

- A full charge enables over 200 times exposure at full power and over 400 times at average power.
- Normal exposure is available even when the power is lower than 20%.




Detachable Battery

Convenient Maintenance & Updates

- Exclusive USB update function supports instant system update
- Different prompt codes indicate system status



5 modes of exposure



App



Hand Switch*



Remote Control



Foot Switch



Delay Key

*Hand Switch is Standard, Others are optional accessories.

Specifications

Max Output Power	5 kW
kV Range	40-125kV
MA's Range	0.4-200mAs
mA	10-100mA



**DIAGNOSTIC
IMAGING SYSTEMS**

QUALITY | INTEGRITY | VALUE...privately owned since 1983

www.vetxray.com



Contact us:

605.341.2433
scott@vetxray.com
2325 E. Saint Charles St.
Rapid City, SD 57703